

ETP electron multipliers

Introduction

Operation

The diagram illustrates the electrical connections for a detector system. It features three main functional blocks connected in series: an 'ion to electron conversion' block, a 'control section', and an 'amplifying section'. The signal path begins with an 'ion input' entering the conversion block, which then passes through the control section and finally the amplifying section, resulting in a 'signal output'. The system is powered by two distinct high-voltage sources. A 'fixed supply ~ 550V' provides a constant potential, with terminals labeled V_1 (connected to the conversion block), V_c (connected to the control section), and V_5 (connected to the amplifying section). A 'variable HV supply for gain control' provides an adjustable potential, labeled V_G , which is connected to the 'mounting surface' of the amplifying section. All components are referenced to a common ground, represented by ground symbols at the bottom of the diagram.

[illegible]

Performance

Conclusions

Figure 10 is a line graph showing the relationship between Measured Counts (normalised at 2850V across amplifying section) and HV across amplifying section of multiplier. The x-axis represents HV across amplifying section of multiplier, ranging from 2200 to 3600 V. The y-axis represents Measured Counts (normalised at 2850V across amplifying section), ranging from 0 to 3. The graph displays four curves corresponding to different attenuation factors:

- Unattenuated, ions = 6.3e4 per sec (Blue diamonds)
- attenuation factor = 90, ions = 6.3e6 per sec (Pink squares)
- attenuation = 8300, ions = 6.3e7 per sec (Red triangles)
- attenuation = 24000, ions = 6.3e8 per sec (Green squares)

The curves show that the measured counts increase with increasing HV across the amplifying section of the multiplier. The unattenuated curve (blue diamonds) shows the highest counts, while the attenuated curves (pink squares, red triangles, green squares) show lower counts, indicating the effect of attenuation on the measured signal.

Figure 10 is a log-linear plot showing Signal Attenuation versus Control Voltage V_c (as percentage of V_r). The y-axis is logarithmic, ranging from $1.0E+05$ to $1.0E+00$. The x-axis is linear, ranging from 0% to 90%. A dashed red line represents the 'No attenuation' case, and a solid red line with markers shows the actual signal attenuation, which increases exponentially with control voltage.

Control Voltage V_c (as percentage of V_r)	Signal Attenuation
0%	$1.0E+05$
10%	$1.0E+04$
25%	$1.0E+03$
45%	$1.0E+02$
65%	$1.0E+01$
85%	$1.0E+00$

ETP electron multipliers

©Copyright 2013 ETP Electron Multipliers. All rights reserved. TP-0064-A RevB