


Installation Notes


Model 14516 ACTIVE FILM Multiplier

ACTIVE FILM Multiplier model 14516 is a direct replacement multiplier for the Hewlett-Packard 5972 MSD (Mass Selective Detector).

NOTE: First Time installation of ACTIVE FILM Multiplier model 14516 may require a connection lead kit (model number LS515). It is included with this multiplier. These leads are reusable and can be used with subsequent 14516 installations.



Bottom view of AF516 showing connection pin positions.



Connection diagram showing leads from the multiplier to the HP5972 flange.